

**MEDIPHOT – LASERFILMS**

Availability of Darkroom (NIF) and Daylight Load Packaging (DLP) for the various cameras.

Brand/Type of Lasercamera	MEDIPHOT Type of Laserfilm		
	X-PHN/RP	X-PLD/RP	X-PIR/RP
Wave Length	633nm	670nm	780-820nm
Packing Specs to be stated when ordering.			
AGFA Matrix Compact L	NIF/DLP-A		
AGFA LR 3300 S	NIF/DLP-A		
AGFA LR 3300 P	NIF/DLP-A		
AGFA 5200 S	NIF/DLP-A		
AGFA 5200 P	NIF/DLP-A		
Sterling LP 300	NIF		
Sterling LP 400	NIF		
Nishimoto ELK EL 2000			NIF
Konica LI-10A	NIF		
Konica LI-21 S4	NIF		
Konica LI-21 D8	NIF		
Kodak Ektascan 100	NIF		
Kodak Ektascan 100 XLP	NIF		
Kodak Ektascan 2180		NIF/DLP-K2	
Kodak Ektascan 1120			NIF/DLP-K1
Kodak Ektascan XLP 400			NIF/DLP-K1
Imation/3M P831			NIF/DLP-M1
Imation/3M M952			NIF/DLP-M1
Imation/3M M959			NIF/DLP-M1
Imation/3M 969 HQ/HS			NIF/DLP-M2
Fuji FL-IM			NIF/DLP-F1

**ART COLENTA GROUP**

**COLENTA Labortechnik  
Ges.m.b.H. & Co. KG**  
Neunkirchner Strasse 117  
A-2700 WIENER NEUSTADT  
Austria  
Tel.: +43/2622/28 311-0  
Fax: +43/2622/28 311-7  
e-mail: [office@colenta.at](mailto:office@colenta.at)

**COLENTA Labortechnik  
Ges.m.b.H. & Co. KG**  
Ratzinger Weg 9  
D-94065 WALDKIRCHEN  
Germany  
Tel.: +49/8581/207-0  
Fax: +49/8581/207-34  
e-mail: [colenta@t-online.de](mailto:colenta@t-online.de)

**ACI - ART COLENTA  
International Ltd.**  
Forge Close, Little End Road,  
Eaton Socon,  
CAMBRIDGESHIRE PE 19 3JH,UK  
Tel.: +44/1480/21 11 18  
Fax: +44/1480/21 52 36  
eFax: +44/870/132 1352  
e-mail: [aci@patflaherty.com](mailto:aci@patflaherty.com)

More info regarding ART COLENTA Medical Product Line  
see <http://www.colenta.at> and <http://www.colenta.de>