

MEDIPHOT

Comparison of X-Ray Film Processors regarding capacity:

1. Films per hour of size 35/43 – 35/35 – 24/30 – 13/18
2. Films per hour of average size approx. 24/30
3. Intake speeds at defined cycle time
4. Waiting time for operator till next film can be fed into processor
5. Total access-time to films to be processed of various sizes

	Film Sizes														
	Min. cycle min	Speed cm/min	Max. width cm	35 x 43 cm			35 x 35 cm			24 x 30 cm			13 x 18 cm		
				Films/hour	Access. Time for proc. Film min	Waiting time for operator till next film can be fed min	Films/hour	Access. Time for proc. Film min	Waiting time for operator till next film can be fed min	Films/hour	Access. Time for proc. Film min	Waiting time for operator till next film can be fed min	Films/hour	Access. Time for proc. Film min	Waiting time for operator till next film can be fed min
Small units offered at low prices	2	40	35	55	3,08	1,08	69	2,88	0,88	100	2,6	0,6	267	2,45	0,45
MEDIPHOT 937	1,5	100	35	140	1,93	0,43	171	1,85	0,35	250	1,74	0,24	667	1,68	0,18
MEDIPHOT 943	1,5	100	43	171	1,85	0,35	171	1,85	0,35	250	1,74	0,24	923	1,63	0,13
MEDIPHOT 902S	1,5	125	43	214	1,78	0,28	214	1,78	0,28	312	1,69	0,19	1.154	1,60	0,10
MEDIPHOT 903	1,5	175	43	300	1,70	0,20	300	1,70	0,20	437	1,64	0,14	1.615	1,57	0,07
										("average")					

COLENTA LABORTECHNIK
 Ges.mbH. & Co. KG
 Neunkirchner Strasse 117
 A-2700 Wiener Neustadt
 Tel.: +43-2622-28311-0
 Fax: +43-2622-28311-7
 E-mail: office@colenta.at
<http://www.colenta.at>

ACI-COLENTA
 INTERNATIONAL Div.
 Capix House, 29 Levellers Lane,
 St. Neots. Cambs. PE19 2JL – UK
 Tel. +44-1480-211118
 Fax: +44-1480-215236
 EMail: +44/870/132 1352
 E-mail: aci@patflaherty.com