



INTRAORAL DENTAL FILM PROCESSOR PERIOMAT

for use in daylight



After approx. 400 films (or 3 weeks of use) the film transport rack and the chemistry tanks can be easily cleaned/rinsed with warm water.

with Daylight Hood

with temp. controlled bath heating and film dryer

Suitable film sizes:

without adapter	3 x 4 cm
with adapter	5,7 x 7,6 cm
	4 x 5 cm
	2,7 x 5,4 cm
	2 x 3,5 cm
	2,4 x 4 cm

Adapters are supplied with the processor.

Tech. specs.:

Voltage	230V
Frequence	50Hz
Total output	400W
Current	1,8A
Duty rating	100%
Process time	5min
Bath temp.	24°C
Weight	11,5kw
Size (LxWxH)	63x25x40cm

COLENTA LABORTECHNIK
Ges.m.b.H. & Co. KG
Neunkirchner Strasse 117
A-2700 Wiener Neustadt
Tel.: +43-2622-28311-0
Fax: +43-2622-28311-7
E-mail: Colenta@nextra.at
<http://www.art-colenta.com>

COLENTA LABORTECHNIK
Ges.m.b.H. & Co. KG
Ratzinger Weg 9
D-94065 ,Waldkirchen
Tel.: +49-8581-207-0
Fax: +49-8581-207-34
E-mail: Colenta@t-online.de
<http://www.colenta.de>

ACI-COLENTA
INTERNATIONAL Div.
Forge Close, Little End Rd,
Eaton Socon, Cambs. UK
Tel. +44-1480-211118
Fax: +44-1480-215236
EFax: +44/870/132 1352
E-mail: aci@patflaherty.com