

DX Easy Imaging Standard CR/DR

The DX Easy Imaging Standard DR is a space-saving and reliable all digital X-ray system. The optimized workflow of this system reduces the number of steps in achieving brilliant images for convenient diagnosis.

Highlights

- Efficient and fast workflow
- Flexible adjustable Z-Arm design
- Functional, space-saving design
- Perfectly suitable for images of skeleton and thorax
- Optional available with detector in cassette size
- Modern touch screen function
- Installation in one day
- Also available as Colenta system solution
- Powered by Software DX Easy Imaging



Technical Specifications



Universal stand	<p>Height adjustment: 27 – 160 cm, electromagnetic brakes</p> <p>Cross arm: rotatable, inclined exposures possible, electromagnetic brakes</p> <p>SID: 120 – 200 cm</p> <p>Bucky unit: ± 90° rotatable</p> <p>Weight: approx. 280 kg (617 lbs), incl. tube, collimator, bucky and base plate</p>
Generator	CMP 200 DR, 50 kW, max. 630 mA, max. 150 kV, max. 200 kHz, AEC, hand switch, 3 phases/N/PE 400 V, incl. console
Tube	Toshiba E7252X, focal spots 0.6 / 1.2 mm, 300 KHu, max. 150 kV, incl. high voltage cable, weight approx. 18 kg (39.7 lbs)
Collimator	Manually turnable, LED lamp, integrated timer, laser alignment
DAP	Dose area product measurement device, integrated display, automatic integration into X-AQS X-ray journal
Measuring Chamber AEC	3-feld carbon measuring chamber with cable
Optional	Tube: RAD-14, RAD-21 Other generators on request
Warranty	12-month warranty Extended warranty on request

COLENTA Labortechnik GmbH & Co KG

Neunkirchner Strasse 117
A-2700 Wiener Neustadt
Austria
Phone: +43-2622-28311-0
E-Mail: office@colenta.at