



Colenta

Medical Product Portfolio


to cover all the needs of Radiologists / X-Ray Departments in Hospitals



**Automatic X-Ray
Film Processors**



**X-Ray Films
Cassettes / Screens**



**CR / DR Digital X-Ray
Retrofit Kit (WiFi)**



**Manual X-Ray Film
Processing Equipment**



**X-Ray Film
Viewers**



**X-Ray
Protection Apparel**

Automatic X-Ray Film Processors * Processing Chemistry *****

X-Ray Film Intensifying Screens * X-Ray Film Cassettes *****

X-Ray Protective Apparel * X-Ray Darkroom Accessories *****

X-Ray Film Viewers * Wireless DR X-Ray Systems *****

Dry Imaging Printer and Films

For further details please visit our homepage www.colenta.at !

COLENTA Labortechnik GmbH & Co. KG

A-2700 Wiener Neustadt Neunkirchner Str. 117 Austria

Tel.: +43(0)2622-28311-0 E-Mail: office@colenta.at

Colenta[®]

MEDIPHOT X-Ray product range

The **Colenta Group** has a world-wide reputation as the Professional Processor Manufacturer, producing equipment and systems for the developing of sensitised film and papers. Our **Mediphot x-ray products** are designed for the automatic developing of x-ray images and include a comprehensive range of x-ray film, chemistry consumables and support items to meet a varied range of x-ray department requirements and applications.

Mediphot products include:

- ✓ a comprehensive range of automatic x-ray film processors (including dental)
- ✓ competitively priced analogue x-ray films, green and blue sensitized
- ✓ including mammography film
- ✓ processing and cleaning chemistry
- ✓ darkroom accessories
- ✓ x-ray protection apparel
- ✓ x-ray cassettes
- ✓ x-ray intensifying screens
- ✓ ID camera systems
- ✓ x-ray film viewers
- ✓ darkroom equipment
- ✓ Retrofit Kits for digital conversion

We want to draw your attention to our DR systems we have launched recently:

The DR System we are offering has an AED (Automatic Exposure Detection), it does not need any interlink to the x-ray generator anymore! You take the DR panel as you take now the x-ray film cassette, you insert the DR panel either in the x-ray table, or in the stand or whatsoever, and the assistant makes the x-ray shoot as he is used to do it with the film. The image is immediately transferred wireless to the computer, and there, our acquisition software creates the image for the Doctor! Our software is well appreciated, as it gives a lot of assistance to the x-ray staff in the Hospital. Our DR panel represents the most modern, latest development in DR technology, and is available with our well-appreciated software!

For further details please visit our homepage www.colenta.at !