





## ANTIX MOBILE PROTECTIONS

	<p>Mobile protection on rubbered wheels with brakes. Stove enamelled steel structure. Antix panel in wood with interposed lead Pb 2 mm. Total coated in plastic laminated. Complete with lead glass and plate pocket.</p> <p><b>MOBILE PROTECTION 1 FLAP</b> Protection 80x190h cm with lead glass 24x30 cm – weight 73 kg Protection 80x190h cm with lead glass 30x40 cm – weight 75 kg</p> <p><b>MOBILE PROTECTION 3 FLAPS</b> Protection 30+80+30x190h cm with lead glass 24x30 cm – weight 105 kg Protection 30+80+30x190h cm with lead glass 30x40 cm – weight 107 kg</p> <p>According to CEE rules applicable to the subject. CE</p>
	<p><b>PROTECTION FOR ANGIOGRAPHY</b></p> <p>Stove enamelled steel structure, on rubbered turning wheels with brakes. Constituted by 2 panels in antix wood Pb 2 mm, one of which zips and complete of lead glass Pb 2 mm 30x40 cm and antix apron Pb 0.50 size 50x30 cm. Panels dimension 55x180h cm + 50x65h cm.</p> <p>According to CEE rules applicable to the subject. CE</p>
	<p><b>ANTIX SCREENS WITH LEAD GLASS</b></p> <p>Constituted by antix panel in the inferior part and by lead glass in upper part for all width. Screening Pb 2 mm. Perimetric aluminium structure on turning rubbered wheels with brakes.</p> <p>Model 87x190h cm with lead glass cm 80x60h</p> <p>Model 107x190h cm with lead glass 100x60h cm</p> <p>According to CEE rules applicable to the subject. CE</p>
	<p><b>MOBILE SCREEN WITH LEAD GLASS</b></p> <p>Heavy stove steel structure. Upper part with lead glass equiv. Pb 2 mm size 50x40 cm and antix apron with rubbered strips equiv. Pb 0.25 mm size 50x25 cm. Lead glass and antix apron adjustable in height. Inferior part with antix panel Pb 1 mm size 100x75h cm.</p> <p>According to CEE rules applicable to the subject. CE</p>

**COLENTA Labortechnik  
GmbH & Co. KG**  
Neunkirchner Strasse 117  
A-2700 Wiener Neustadt  
Austria  
Tel: +43/2622/28 311-0  
Fax: +43/2622/28 311-7  
e-mail: [office@colenta.at](mailto:office@colenta.at)

**ACI – ART COLENTA  
International Ltd.**  
Forge Close, Little End Road,  
Eaton Socon,  
CAMBRIDGESHIRE PE 19 3JH, UK  
Tel: +44/1480/21 11 18  
Fax: +44/1480/21 52 36  
e-mail: [aci@patflaherty.com](mailto:aci@patflaherty.com)

More info regarding COLENTA Medical Product Line  
see <http://www.colenta.at>